

DISPOSABLE
HEALTHCARE
PRODUCTS



Keeping pace with innovation

EXEL^{INT}
ANIMAL HEALTH DIVISION



CONTENTS

Insulin Syringes	1
U- 40 Insulin Syringes	2
Tuberculin Syringes	3
3cc Syringes	4
5-6cc Syringes	5
10-12cc Syringes	5
Curved Tip Syringes	5
20-25cc Syringes	6
30-35cc Syringes	6
50-60cc Syringes	6
Culinary Syringe	6
Aluminum Hub Needles : Hard Pack	7
Aluminum Hub Blunt Needles	7
Hypodermic Needles	8
Drenching Syringes	9
Bulb Syringes	9
Multi-Sample Needles	10
Multi-Sample Luer Adapters	10
Multi-Sample Luer Holders	10
Injection Plugs	10
Vaculet Blood Collection Sets	11
Scalp Vein (Butterfly) Sets	11
IV Catheters	12
Spinal Needles	13
IV Administration Sets	14
Surgical Blades	15
Disposable Scalpels	15
Bulk Non-Sterile Syringes	16
Product Specifications	17

Comfort point™ INSULIN SYRINGES

0.3cc LO-DOSE INSULIN SYRINGE

with permanently
attached needle

Catalog No. Description

26018	U-100 29G x 1/2"
26014	U-100 30G x 5/16"

10/bag,100/box,500/case

1/2cc LO-DOSE INSULIN SYRINGE

with permanently
attached needle

Catalog No. Description

26026	U-100 28G x 1/2"
26028	U-100 29G x 1/2"
26015	U-100 30G x 5/16"

10/bag,100/box,500/case

1cc INSULIN SYRINGE

with permanently
attached needle

Catalog No. Description

26027	U-100 28G x 1/2"
26029	U-100 29G x 1/2"
26016	U-100 30G x 5/16"

10/bag,100/box,500/case

- This product does not contain natural rubber latex.
- Sterile
- Non-toxic
- Non pyrogenic
- Permanently attached needle
- Reduced dead space for reduced Insulin waste
- Smooth plunger action for greater filling accuracy
- Auto ground needle tip for greater comfort



Comfort point™

U - 40 INSULIN SYRINGES



- This product does not contain natural rubber latex.
- Sterile
- Non-toxic
- Non pyrogenic
- Permanently attached needle
- Reduced dead space for reduced Insulin waste
- Smooth plunger action for greater filling accuracy
- Auto ground needle tip for greater comfort

EXEL ANIMAL HEALTH DIVISION

EXEL U-40 *Comfort point* Sterile Insulin Syringes

Comfort point a name recognized for quality at affordable pricing.
 Ask your EXEL distributor for EXEL's U-40 Comfort Point Insulin Syringes ...
 If it's not Comfort Point U-40 it's not EXEL. ... and start saving now.

EXEL is your proven partner in animal healthcare.

Catalog No.	Description
26036	U-40 1/2cc 29G x 1/2"
26037	U-40 1cc 29G x 1/2"

10/bag, 100/box, 500/case

with permanently attached needle

TUBERCULIN SYRINGES

1cc TUBERCULIN SYRINGE*

with detachable needle, Luer-slip

Catalog No. Description

26040*	1cc 27Gx1/2"
26042*	1cc 25Gx1"
26043*	1cc 26Gx1/2"
26044*	1cc 25Gx5/8"

100/box, 1,000/case



*Low Dead Space Plunger



Luer Slip

Luer Lock

1cc TUBERCULIN SYRINGE ONLY

Catalog No. Description

26048*	1cc Slip Tip with cap
26049*	1cc Lock Tip with cap
26050	1cc Lock Tip

100/box, 1,000/case

1cc ZERO DEAD SPACE TUBERCULIN SYRINGE

with permanently attached needle

Catalog No. Description

26046	1cc 25Gx5/8"
-------	--------------

10/bag, 100/box, 500/case

- This product does not contain natural rubber latex.
- Sterile
- Non-toxic
- Non pyrogenic
- Reduced dead space for reduced medication waste
- Smooth plunger action for greater filling accuracy
- Auto ground needle tip for greater comfort



3CC SYRINGES



3cc SYRINGE only

Catalog No. Description

26200	Luer Lock with cap
26201	Luer Slip with cap
100/box, 1,000/case	

- This product does not contain natural rubber latex.
- Sterile
- Non-toxic
- Pyrogen free
- Low Dead Space
- With Cap
- Positive plunger stop
- Crystal clear barrels
- Easy to read large, bold print



LUER -LOCK TIP

3cc Syringe/Needle Combination

Catalog No.	Description	Hub Color
26100	25G x 5/8"	Orange
26101	23G x 1"	Blue
26102	22G x 1"	Black
26103	22G x 1 1/4"	Black
26104	22G x 1 1/2"	Black
26105	21G x 1 "	Green
26107	21G x 1 1/2"	Green
26108	20G x 1 "	Yellow
26109	20G x 1 1/2"	Yellow
26110	18G x 1 1/2"	Pink
26111	25G x 1 "	Orange
26112	25G x 1 1/2"	Orange
26115	22G x 3/4"	Black
26117	23G x 1 1/2"	blue
26127	21G x 1 1/4"	Green
26129	27G x 1 1/4"	Medium Grey

100/box, 1,000/case

LUER -SLIP TIP

3cc Syringe/Needle Combination

Catalog No.	Description	Hub Color
26601	22G x 3/4"	Black
26602	25G x 5/8"	Orange
26603	23G x 1"	Blue
26605	22G x 1 "	Black
26608	21G x 1 "	Green
26611	20G x 1 "	Yellow

100/box, 1,000/case

5-6CC SYRINGES 10-12CC SYRINGES CURVED TIP SYRINGES

5-6cc Syringe/ Needle Combination LUER -LOCK TIP

Catalog No.	Description	Hub Color
26210	22G x 1"	Black
26211	22G x 1 ½"	Black
26212	21G x 1"	Green
26213	21G x 1 ½"	Green
26214	20G x 1"	Yellow
26218	21G x 1 1/4"	Green
26219	22G x 1 1/4 "	Black

100/box, 800/case

5-6cc SYRINGE only

Catalog No.	Description
26230	Luer Lock with cap
26231	Luer Slip with cap

100/box, 800/case



- This product does not contain natural rubber latex.
- Sterile
- Non-toxic
- Pyrogen free
- Positive plunger stop to prevent spillage when drawing
- Smooth plunger motion

10-12cc Syringe/ Needle Combination LUER -LOCK TIP

Catalog No.	Description	Hub Color
26250	22G x 1"	Black
26251	22G x 1 ½"	Black
26252	21G x 1"	Green
26253	21G x 1 ½"	Green
26254	20G x 1"	Yellow
26255	20G x 1 ½"	Yellow
26257	18G x 1 "	Pink

100/box, 800/case

10-12cc SYRINGE only

Catalog No.	Description
26265	Luer Lock with cap
26266	Luer Slip with cap

100/box, 800/case

CURVED TIP SYRINGE

Catalog No.	Description
26267	Curved Tip

50/box, 500/case



20-25CC SYRINGES 30-35CC SYRINGES 50-60CC SYRINGES CULINARY SYRINGE

without needle

20-25cc SYRINGE only

Catalog No. Description

26280 Luer Lock with cap
26281 Luer Slip Ecc with cap
50/box, 500/case

30-35cc SYRINGE only

Catalog No. Description

26290 Luer Lock with cap
26291 Luer Slip Ecc with cap
26292 Cath Tip Ecc
50/box, 250/case

50-60 SYRINGE only

Catalog No. Description

26300 Luer Lock with cap
26301 Luer Slip Ecc with cap
26304 Catheter Tip with cap
25/box, 150/case

50-60cc SYRINGE only

Bulk, non sterile, lab clean

Catalog No. Description



Bn26302 Catheter Tip with cap
400/case

CULINARY SYRINGE

Catalog No. Description

26140 Culinary Syringe
with Culinary needle
100/bag

- This product does not contain natural rubber latex.
- Sterile
- Non-toxic
- Pyrogen free
- Positive plunger stop to prevent spillage when drawing
- Smooth plunger motion
- Polypropylene construction



ALUMINUM HUB NEEDLES

ALUMINUM HUB BLUNT NEEDLES

ALUMINUM HUB HYPODERMIC NEEDLES : HARD PACK

Catalog No.	Description
26450	14G x 1"
26452	14G x 1 1/2"
26454	14G x 2"
26455B	16G x 5/8"
26456	16G x 3/4"
26457	16G x 1"
26458	16G x 1 1/2"
26461	18G x 3/4"
26462	18G x 1"
26463B	18G x 5/8"
26464	18G x 1 1/2"
26467	20G x 1"
26468	20G x 1 1/2"

100/box,1,000/case



- Designed for veterinary use
- Sterile, Non toxic
- Pyrogen free
- Autoclavable
- Color-coded for easy identification
- Each needle is individually packaged and sealed with a tamper evident heat stake.

ALUMINUM HUB BLUNT NEEDLES

Catalog No.	Description
26471	15G x 1 1/2"
26472	17G x 1 1/2"
26473	18G x 1"
26474	19G x 1 1/2"
26475	20G x 1 1/2"

25/box,100/case



- This product does not contain natural rubber latex.
- Sterile, single-use, individually packaged blunt needle
- Stainless steel laser welded, non-lubricated cannula
- Aluminium hub, luer lock design
- Hub crimp cannula attachment
- Color-coded for easy identification
- 2-piece polypropylene needle cap and sheath



HYPODERMIC NEEDLES



HYPODERMIC NEEDLES REGULAR BEVEL

Catalog No.	Description	Hub Color
26400	27G x 1/2"	Medium Grey
26402	26G x 1/2"	Brown
26403	25G x 5/8"	Orange
26405	25G x 1"	Orange
26406	25G x 1 1/2"	Orange
26407	23G x 3/4"	Blue
26408	23G x 1"	Blue
26409	23G x 1 1/2"	Blue
26410	22G x 3/4"	Black
26411	22G x 1"	Black
26412	22G x 1 1/2"	Black
26414	21G x 1"	Green
26416	21G x 1 1/2"	Green
26417	20G x 1"	Yellow
26418	20G x 1 1/2"	Yellow
26419	18G x 1"	Pink
26420	18G x 1 1/2"	Pink
26424	19G x 1"	Cream
26426	27G x 1 1/2"	Medium Grey
26427	27G x 1 1/4"	Medium Grey
26430	20G x 3/4"	Yellow
26431	20G x 1/2"	Yellow
26438	19G x 1 1/2"	Cream



SPECIALITY USE HYPODERMIC NEEDLES

Catalog No.	Description	Hub Color
26437	30G x 1/2"	Yellow
26439	30G x 1"	Yellow
100/box, 2,000/case		
26448	28G x 3/4"	Blue Green
100/box, 1,000/case		
26391	31G x 1/2"	Orange
26392	32G x 1/2"	Green
26392S	32G x 1/2"	Green Soft Pack
26393	33G x 1/2"	Clear
100/box/1000/case		

100/box, 2,000/case

NEEDLE LENGTH

Inches	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"	2"
Millimeters	13mm	16mm	19mm	25mm	32mm	38mm	50mm

- This product does not contain natural rubber latex.
- Sterile, Non-toxic, Pyrogen free
- SUS - 304 stainless steel
- Needles and boxes are color-coded for easy identification
- Each needle is individually packaged, sterilized and ready for use

DRENCHING SYRINGES

BULB SYRINGES

FEEDING TUBES

- This product does not contain natural rubber latex.
- Made of resilient Thermoplastic material.
- Color coded for easy identification
- Raised black printing
- Replacement O-rings included
- Syringes packaged in attractive peg-able packaging



DRENCHING SYRINGES

Catalog No.	Description	Hub Color
25500	10ml Short (for Small Hands)	Yellow
25502	10ml	Blue
25504	20ml	Red
25506	30ml	Green
25510	Feeding Nozzle	

12/box

- Useful in the kidding barn for clearing airways of newborns
- Useful in irrigating ears or small cavities
- Flat bottoms allow syringe to stand up without tipping
- Complete manual control
- Bulk, Non-sterile

BULB SYRINGES

Catalog No.	Description
25531	3 oz

50/Case



MULTI-SAMPLE NEEDLES

MULTI-SAMPLE LUER ADAPTERS

MULTI-SAMPLE LUER HOLDERS

INJECTION PLUG



MULTI-SAMPLE (BLOOD-DRAW) NEEDLES

Catalog No.	Description	Hub Color
26500	18G x 1"	Pink
26501	20G x 1"	Yellow
26502	20G x 1 ½"	Yellow
26503	21G x 1"	Green
26504	21G x 1 ½"	Green
26505	22G x 1"	Black
26506	22G x 1 ½"	Black

100/box, 1000/case



- This product does not contain natural rubber latex.
- Sterile
- Non-toxic
- Multi-sample needles permit several samples to be taken with a single puncture
- Sharp & smooth edges make penetration painless for patients and provides easy connection to rubber stoppers

MULTI-SAMPLE LUER HOLDERS

Catalog No.	Description
26530	5-7.5ml holder

100/bag, 1000/case (Non-Sterile)

INJECTION PLUG

Catalog No.	Description
26539	Injection plug with cap LL


50/box, 200/case



MULTI-SAMPLE LUER ADAPTERS

Catalog No.	Description
26537	Adapter

100 pcs/box, 1000 pcs/case




BLOOD COLLECTION SETS SCALP VEIN SETS



BLOOD COLLECTION SET

Catalog No.	Description	Tube	Hub Color
26764	21G x 3/4"	12"	Green
26766	23G x 3/4"	12"	Blue
26768	25G x 3/4"	12"	Orange

50/box 200/case



SCALP VEIN SET



Holder not included

PSV / BUTTERFLY / SMALL VEIN INFUSION

Catlog No.	Description	Tube	Hub Color
26702	19G x 3/4"	12"	Cream
26703	20G x 3/4"	12"	Yellow
26704	21G x 3/4"	12"	Green
26706	23G x 3/4"	12"	Blue
26708	25G x 3/4"	12"	Orange
26709	27G x 3/4"	12"	Medium Grey

50/box,500/case



- This product does not contain natural rubber latex.
- Individually packaged
- Thin wall stainless steel cannula, siliconized for easy insertion
- Flexible wings hold the needle in position and help to prevent inadvertent pull out
- Color-coded for easier identification

SAFELET IV CATHETERS



FEATURES

- Easy dispenser pack
- Functional design casing permitting one action opening and no rolling
- Color-coded cap allows for easier identification of catheter size
- Translucent catheter hub and flashback chamber allows for easy detection of blood flashback at vein insertion
- Radio-opaque striped catheter
- Precision finished FEP catheter assures stable flow and eliminates catheter tip kink during venipuncture
- Cannula with additional back bevel
- Can be connected to syringe by removing filter cap to expose Luer taper end
- Use of hydrophobic membrane filter eliminates blood leakage
- Close and smooth contact between catheter tip and inner needle enables safe and smooth venipuncture
- Crimped smooth engagement between catheter and hub, without use of adhesive, helps eliminate blood residuals and helps reduce blood coagulation or hemolysis

This product does not contain natural rubber latex.

SPECIFICATIONS


Item No.	Size	Stainless Cannula		Hub Color
		Gauge	Outside Diameter (mm)	
26726	14G x 1 1/4"	14G	1.63	Orange
26729	14G x 2"	14G	1.63	Orange
26730	16G x 1 1/4"	16G	1.2	Medium Grey
26731	16G x 2"	16G	1.2	Medium Grey
26735	18G x 1 1/4"	18G	0.9	Green
26736	18G x 1 1/2"	18G	0.9	Green
26738	18G x 2"	18G	0.9	Green
26741	20G x 1"	20G	0.71	Pink
26742	20G x 1 1/4"	20G	0.71	Pink
26743	20G x 2"	20G	0.71	Pink
26746	22G x 1"	22G	0.56	Blue
26751	24G x 3/4"	24G	0.41	Yellow



50/ box, 200/case

SPINAL NEEDLES

- This product does not contain natural rubber latex.
- Convenient blister packaging
- Siliconized for smooth comfortable penetration
- Transparent plastic hub



'QUINKE' BEVEL SPINAL NEEDLES

Catalog No.	Description	Hub Color
26960	18G (1.20mm) x 3 1/2"	Pink
26961	20G (1.20mm) x 2 3/4"	Yellow
26964	20G (0.90mm) x 3 1/2"	Yellow
26965	22G (0.70mm) x 2 3/4"	Black
26967	22G (0.70mm) x 3 1/2"	Black
26970	25G (0.50mm) x 3 1/2"	Orange

50/ box, 100/case




IV ADMINISTRATION SETS

I.V. ADMINISTRATION SET VENTED/NON VENTED

With Y-injection Site, Vented / Non Vented, With Clamp • 78" Tubing

Priming Volume
23ml/Approx
Inside Diameter
3.0 ± 0.05mm
Outside Diameter
4.1 + 0.10 / -0.05



Catalog No. Description

29081 15 Drop/ml Luer Slip
29083 60 Drop/ml Luer Slip

50 pcs/case



Catalog No. Description

29082 15 Drop/ml Luer Lock
29085 60 Drop/ml Luer Lock

50 pcs/case



Catalog No. Description

29086 15 Drop/ml Lock/Slip
27086 15 Drop/ml Lock/Slip with Needle-Free Injection Site

50 pcs/case

I.V. ADMINISTRATION SET VENTED/NON VENTED 105"

Priming Volume
28ml/Approx
Inside Diameter
3.0 ± 0.05mm
Outside Diameter
4.1 + 0.10 / -0.05




Catalog No. Description

29088 15 Drop/ml Lock/Slip, Vented/Non Vented, 2 Injection Sites,
Roller Clamp, 2 Occlusion Clamps

50 pcs/case

- This product does not contain natural rubber latex.
- Sterile
- Non toxic
- Non pyrogenic fluid path
- Vented/ Non Vented
- Roller clamp
- Occlusion clamp
- Pouch packaging

SURGICAL BLADES DISPOSABLE SCALPELS



- This product does not contain natural rubber latex.
- Made from premium quality Swiss Carbon Steel and British Stainless Steel
- Computer controlled grinding, sharpening and polishing, ensures consistently smooth incision properties
- Sterilized by GAMMA Radiation 2.5m - RADS
- A Vapor Corrosion Inhibitor (VCI) paper lining, leakage testing, and an aluminum foil peel pack help insure quality and integrity when ready for use
- Each production lot is sampled and conforms to ISO 2859, inspection level (S4) AQL 1.0 & 2.5
- Blade rigidity is tested to eliminate blade breakage
- Blades conform to ISO: 7740 (1995) & BS: 2982 (1986)
- Blades are ultrasonically cleaned prior to packing

STAINLESS STEEL SURGICAL BLADES



Catalog No.	Size
29500	10
29502	11
29504	12
29508	15
29515	22
100/box, 1,000/case	

CARBON STEEL SURGICAL BLADES

Catalog No.	Size
29526	11
29532	15
100/box, 1,000/case	

STAINLESS STEEL DISPOSABLE SCALPELS

Catalog No.	Size
29550	10
29552	11
29556	15
29559	22
10/box, 100/case	



BULK NON-STERILE SYRINGES



This product does not contain natural rubber latex.

1cc BULK NS SYRINGE

Catalog No.	Description	Qty. Per Case
BN 26048	Luer Slip	1000
BN 26050	TB Luer Lock	4000

3cc BULK NS SYRINGE

Catalog No.	Description	Qty. Per Case
BN 26200	Luer Lock	1500
BN 26201	Luer Slip	1500

5-6cc BULK NS SYRINGE

Catalog No.	Description	Qty. Per Case
BN 26230	Luer Lock	2000
BN 26231	Luer Slip	2000

10-12cc BULK NS SYRINGE

Catalog No.	Description	Qty. Per Case
BN 26265	Luer Lock	1500
BN 26266	Luer Slip	1500

20-25cc BULK NS SYRINGE


Catalog No.	Description	Qty. Per Case
BN 26280	Luer Lock	1000
BN 26281	Luer Slip	1000

30-35cc BULK NS SYRINGE

Catalog No.	Description	Qty. Per Case
BN 26290	Luer Lock	400
BN 26291	Luer Slip	400
BN 26292	Cath Tip	400

50-60cc BULK NS SYRINGE

Catalog No.	Description	Qty. Per Case
BN 26300	Luer Lock	400
BN 26302	Cath Tip	400



Luer Slip Eccentric

Catheter Tip



Luer Slip

Luer Lock



* Low Dead Space Plunger



PRODUCT SPECIFICATIONS

PRODUCT	QUANTITY WEIGHT			INCHES			MILLIMETER		
	PCS	LBS	KGS	LENGTH	WIDTH	HEIGHT	LENGTH	WIDTH	HEIGHT
1/2-1cc Insulin	500	4.82	2.19	17	7	6	432	178	152
1cc Insulin	500	5.36	2.43	17	7	6	432	178	152
0.3cc Insulin	500	4.50	2.04	17	7	6	432	178	152
1cc TB 0 Dead - Space	500	5.38	2.44	17	7	7	432	178	178
1cc TB w/needle	1000	14.50	6.58	23	16	10	584	406	254
1cc TB Only	1000	11.50	5.22	19	13	10	483	330	254
3cc w/Needle	1000	16.45	7.46	28	16	9	711	406	229
3cc only	1000	14.30	6.49	28	12	11	711	305	279
5-6cc w/needle	800	17.60	7.98	27	16	9	686	406	229
5-6cc only	800	15.20	6.89	18	15	12	457	381	305
10-12cc w/needle	800	26.40	11.97	30	18	12	762	457	305
10-12cc only	800	21.80	9.89	30	14	12	762	356	305
20-25cc only	500	26.65	12.09	28	16	15	711	406	381
30-35cc only	250	17.65	8.01	31	15	10	787	381	254
50-60cc only	150	15.60	7.08	18	17	13	457	432	330
Scalp Vein Set Butterfly	500	7.00	3.18	18	12	10	457	305	254
Blood Collection Set	200	3.30	1.50	12	10	8	305	254	203
IV Catheter	1000	3.30	1.50	11	6	6	279	152	152
Hypodermic Needles	2000	6.55	2.97	17	14	5	432	356	127
Multi-Sample Needle	1000	4.60	2.09	16	7	4	406	178	102
Luer Adaptor	1000	5.37	2.44	11	7	7	279	178	178
Luer Holder	1000	6.30	2.86	15	13	12	381	330	305
Spinal Needle	1000	2.25	1.02	12	9	5	305	229	127
IV Administration Set	50	4.31	1.95	15	9	8	381	229	203
Surgical Blades	1000	2.74	1.24	8	7	3	203	178	76
Disposable Scalpels	1000	2.48	1.12	10	9	4	254	229	102
Injection Plugs w/cap	200	1.56	0.71	12	8	3	305	203	76
Drenching Syringes	12	1.54	0.70	2	12	8	51	305	203
Aluminum Hub Needles	1000	5.82	2.64	10	9	5	254	229	127
Aluminum Hub Blunt Needles	100	1.00	0.45	8	3	2	203	76	51
Bulb Syringes	50	6.30	2.86	14	14	5	356	356	127

Specifications vary slightly between sizes.
Packaging subject to change.
Please contact us for information on products not listed here.



EXELINT INTERNATIONAL, CO.

Corporate Headquarters

2500 Santa Fe Avenue

Redondo Beach

California 90278 USA

Phone : (310) 649-0707

Fax : (310) 649-1178

E-mail : info@exelmed.com

Sales & Marketing

1510 Dr. MLK Jr. Street North

Saint Petersburg,

Florida 33704 USA

Phone : (800) 940-3935

Fax : (800) 308-5048

FOR INFO PLEASE CALL

1-800-940-EXEL(3935)

info@exelint.com

www.exelint.com